

Biliary Stent

Used to drain bile from obstructed biliary duct.

Technical information:

- Diameter – 7 fr./10 fr.
- Length – 5cm/15cm.
- Material – Radio-opaque PTFE or PU.
- Variants
 - Straight – S
 - Double pig tail – DP
 - Amsterdam – AMST
 - Standard - STD
- **Sterile for Single use only.**
- Wire guide compatibility – 0.038”.



Ordering information – Mention product code, length and quantity required

Reference	Product Name	Variant	Material
SBS – TEF – 7 - S	Biliary Stent 7 Fr.	Straight	PTFE
SBS – TEF – 10 - S	Biliary Stent 10 Fr.	Straight	PTFE
SBS – PU – 7 - S	Biliary Stent 7 Fr.	Straight	PU
SBS – PU – 10 - S	Biliary Stent 10 Fr.	Straight	PU
SBS – PU – 7 - DP	Biliary Stent 7 Fr.	Double Pig Tail	PU
SBS – PU – 10 - DP	Biliary Stent 10 Fr.	Double Pig Tail	PU
SBS – TEF – 7 - STD	Biliary Stent 7 Fr.	Standard	PTFE
SBS – TEF – 10 - STD	Biliary Stent 10 Fr.	Standard	PTFE
SBS – PU – 7 - STD	Biliary Stent 7 Fr.	Standard	PU
SBS – PU – 10 - STD	Biliary Stent 10 Fr.	Standard	PU
SBS – TEF – 7 - AMST	Biliary Stent 7 Fr.	Amsterdam	PTFE
SBS – TEF – 10 - AMST	Biliary Stent 10 Fr.	Amsterdam	PTFE
SBS – PU – 7 - AMST	Biliary Stent 7 Fr.	Amsterdam	PU
SBS – PU – 10 - AMST	Biliary Stent 10 Fr.	Amsterdam	PU

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181