


Aspiration Needle

Used for Endoscopic aspiration into lungs and bronchile.

Technical information:

- Needle length – 15cm..
- Working Length – 110cm.
- Sharp and short bevel. 14,16 and 19 swg.
- Compatible with 2.0/2.8mm instrument channel.
- **Sterile for Single use only.**



Ordering information – Mention product code and quantity required

Reference	Product Name	Needle Gauge	Working Length	Required Channel
AN – 14	Aspiration Needle	14G	110cm	2.0/2.8mm
AN – 16	Aspiration Needle	16G	110cm	2.0/2.8mm
AN – 19	Aspiration Needle	19G	110cm	2.0/2.8mm

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises