

Biopsy Forcep (Bronchoscopy)

Used endoscopically to obtain Mucosal tissue samples from bronchile for microscopic-histopathological examination.

Technical information:

- Type
 - Without spike
- Length – 110cm.
- Compatible with 2.0/2.8mm instrument channel.
- Individually packed.
- **Sterile for Single use only.**



Ordering information – Mention product code and quantity required

Reference	Product Name	Type	Required Channel	Working Length
S – BF – 2	Biopsy Forcep	No Spike	2.0mm	110cm
S – BF – 2.8	Biopsy Forcep	No Spike	2.8mm	110cm

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises