


## Cleaning Brush (Bronchoscopy)

Technical information:

- Type – Both end bristles.
- Shaft – Plastic coated wire.
- Semi soft and channel friendly bristles.
- Compatible with 2.0/2.8mm instrument channel.
- Length – 110cm.
- Atrumatic tip.
- Individually packed.
- **Non-Sterile** for **Single use only**.



Ordering information – Mention product code and quantity required

Reference	Product Name	Required Channel	Working Length
S – CB – 2	Cleaning Brush (Bronchoscopy)	2.0mm	110cm
S – CB – 2.8	Cleaning Brush (Bronchoscopy)	2.8mm	110cm

Contact:-

[www.makendoscopy.com](http://www.makendoscopy.com)

[mak\\_enter@yahoo.com](mailto:mak_enter@yahoo.com)

[mak\\_enter62@rediffmail.com](mailto:mak_enter62@rediffmail.com)

Phone - +91 22 25289181

MAK Enterprises