


Cytology Brush (Bronchoscopy)

Used to collect cells.

Technical information:

- Length – 110cm.
- Required channel diameter 2.0/2.8mm.
- **Sterile for Single use only.**



Ordering information – Mention product code and quantity required

Reference	Product Name	Required Channel
S – CY - 2	Cleaning Brush (Bronchoscopy)	2.0mm
S – CY – 2.8	Cleaning Brush (Bronchoscopy)	2.8mm

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises