

Foreign Body Holding Forcep 2 Prong, 3Prong (Bronchoscopy)

Used to retrieve foreign body from Esophagus.

Technical information:

- Compatible with 2.0/2.8mm instrument channel.
- Working length 110cm.
- Excellent and secure grip to retrieve foreign body from esophagus.
- 2 Prong and 3 Prong.
- **Non-Sterile** for **Multiple use**.



Ordering information – Mention product code and quantity required

Reference	Product Name	Working Length	Required Channel	Jaws
SFBH - 2	F.B. Holder 2 prong	110cm	2.0mm	2 Prong
SFBH - 3	F.B. Holder 3 prong	110cm	2.0mm	3 Prong
SFBH - 2	F.B. Holder 2 prong	110cm	2.8mm	2 Prong
SFBH - 3	F.B. Holder 3 prong	110cm	2.8mm	3 Prong

Range of Foreign Body Holders – 2 prong, 3 prong, Loop Basket, Retrieval Basket, Snare, Shark Tooth Forcep, Alligator Forcep, Monofilament Memory Basket.

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises