

Cleaning Brush

Used to clean endoscope channel.

Technical information:

- Type – Both end bristles.
- Shaft – Plastic coated wire.
- Semi soft and channel friendly bristles for non invasive cleaning of endoscope channel.
- Compatible with 2.8mm instrument channel.
- Length – 200cm.
- Individually packed.
- **Non-Sterile for Single use only.**



Ordering information – Mention product code and quantity required

Reference	Product Name	Required Channel	Working Length
S – CB – DE	Cleaning Brush Double Ended	2.8mm	200cm
S – CB – DE - BS	Cleaning Brush Double Ended Big-Small	2.8mm	200cm
S – CB – P	Port Cleaning Brush	-	-
S – CB – B	Body Cleaning Brush	-	-

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises