

ERCP Cannula

Used endoscopically to opacify the biliary duct with radiopaque contrast material for fluroscopic view.

Technical information:

- Type – Tapered tip/ Ball tip/ Double lumen.
- Catheter tip diameter – 5.5 Fr./6 Fr.
- Recommended wire guide diameter 0.035”
- Length – 200cm.
- Individually packed.
- **Sterile for Single use only.**



Ordering information – Mention product code and quantity required

Reference	Product Name	Tip	Working Length	Required Channel
S – ERC – T - SL	Single Lumen ERCP Cannula Tappered Tip	5.5Fr	200cm	2.8mm
S – ERC – B - SL	Single Lumen ERCP Cannula Ball Tip	5.5Fr	200cm	2.8mm
S – ERC – T - DL	Double Lumen ERCP Cannula Tappered Tip	6.0Fr	200cm	2.8mm

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises