

# Endo Lithotripter

Used for the purpose of C.B.D Stone Retrieval through the scope and for Lithotripsy to Crush Stones in Biliary Duct.

Technical information:

- Working Length – 210cm.
- Compatible with instrument channel 4.2mm.
- Basket diameter – 35mm.
- Basket shape – hexagonal.
- Flush port for contrast.
- Atraumatic bullet type basket tip for ease of insertion and procedure
- **Non-Sterile** for **Multiple use**.
- Sterile lower assembly supply on request.



Ordering information – Mention product code and quantity required

Reference	Product Name	Sheath	Sheath Length
S – EL - SET	Endo Lithotripter Complete Set with Handle	3.7mm	210cm
S - EL – SHEATH	Endo Lithotripter Flexible Sheath with Basket Filament Assembly	3.7mm	210cm
S – EL - B	Endo Lithotripter Basket only	-	-

Contact:-

[www.makendoscopy.com](http://www.makendoscopy.com)

[mak\\_enter@yahoo.com](mailto:mak_enter@yahoo.com)

[mak\\_enter62@rediffmail.com](mailto:mak_enter62@rediffmail.com)

Phone - +91 22 25289181

MAK Enterprises