


Hydro Neg Guide Wire

Used to gain access through difficult biliary and pancreatic strictures to facilitate exchange of compatible wire guided accessories to the desired biliary or pancreatic ducts.

Technical information:

- Size diameter – 0.032”/0.035”/0.038”.
- Length – 450cm.
- Tip – Straight extra filiform.
- Individually packed.
- **Sterile for Single use only.**



Ordering information – Mention product code and quantity required

Reference	Product Name	Length	Wire
S – HYNeg - 03245	HYDRO NEG GUIDE WIRE	450cm	0.032”
S – HYNeg - 03545	HYDRO NEG GUIDE WIRE	450cm	0.035”
S – HYNeg - 03845	HYDRO NEG GUIDE WIRE	450cm	0.038”

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises