


Manual Lithotripter

Used for C.B.D stone lithotripsy.

Technical information:

- Sheath tip – inner cutting edges to perform lithotripsy. Outer tip rounded atraumatic.
- Secure and firm fitting of S.S with handle to facilitate retrieval – lithotripsy procedure.
- Handle – light weight with unidirectional movement.
- **Non-Sterile for Multiple use.**



Ordering information – Mention product code and quantity required

Reference	Product Name	Sheath	Sheath Length
SML	Manual Lithotripter	4.7mm	95cm

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises