


Retrieval Basket

Used to retrieve foreign bodies and C.B.D stones.

Technical information:

- Working Length – 230cm.
- Compatible with instrument channel 2.8mm/3.2mm.
- Basket Diameter – 30mm/35mm/40mm.
- Multifilament, Hexagonal shaped.
- Atraumatic basket tip for ease of insertion and procedure.
- User friendly handle for Lithotripsy.
- Compatible with Lithotripter.
- **Sterile for Single use only.**



Ordering information – Mention product code and quantity required

Reference	Product Name	Basket	Required Channel
S – RB - 4302	Retrieval Basket	30mm	2.8mm
S – RB - 4352	Retrieval Basket	35mm	2.8mm
S – RB - 4402	Retrieval Basket	40mm	2.8mm

Strand diameter – 0.41mm, Working length – 230cm.

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises