

# Self Expandable Metallic Biliary Stent

Used for Stenosis or Obstructed Biliary tract caused by benign or malignant changes.

Technical information:

- Wire introducer system.
- Radiotherapy compatible.
- Expanded Stent Diameter – 10mm.
- Stent Wire Diameter – 0.16mm.
- Stent Length – 60mm/80mm/100mm.
- Uncovered Stent.
- Stent introducer system diameter – 8Fr.
- Compatible with 0.038inch wire guide and instrument channel size – 3.2mm.
- **Sterile for Single use only.**



Ordering information – Mention product code and quantity required

Reference	Product Name	Stent Length	Stent
S – MBS - 6	Self Expandable Metallic Biliary Stent	60mm	10mm
S – MBS - 8	Self Expandable Metallic Biliary Stent	80mm	10mm
S – MBS - 10	Self Expandable Metallic Biliary Stent	100mm	10mm

Contact:-

[www.makendoscopy.com](http://www.makendoscopy.com)

[mak\\_enter@yahoo.com](mailto:mak_enter@yahoo.com)

[mak\\_enter62@rediffmail.com](mailto:mak_enter62@rediffmail.com)

Phone - +91 22 25289181

MAK Enterprises