


Sphincterotome

Used for Endoscopic Cannulation of Biliary Ductal System and Sphincterotomy.

Technical information:

- Triple Lumen with Precurved tapered tip.
- Working Length – 200cm.
- Cutting wire with protective insulation sheath, unique protective feature for endoscope.
- Unique guide wire port to pass the wire straight without any curve and less friction on guide wire coating.
- Active cord available seperately.
- **Sterile for Single use only.**



Ordering information – Mention product code and quantity required

Reference	Product Name	Tip Length	Working Length	Insulated Wire Length	Cutting Wire	Total Wire Length	Wire Guide Compatibility
SSPH - SN	Sphincterotome Short Nose	6mm	200cm	8mm	S.S.Monofilament	25mm	0.035"
SSPH - LN	Sphincterotome Long Nose	10mm	200cm	8mm	S.S.Monofilament	25mm	0.035"
SSPH - NN	Sphincterotome No Nose	0mm	200cm	0mm	S.S.Monofilament	15mm	0.035"

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181