


Wire Guided Stone Extraction Balloon Catheter

Used for Endoscopic extraction of Biliary Stones.

Technical information:

- Double/Triple lumen feature provides wire guide access to facilitate stone extraction, from desired location during procedure.
- Balloon Length – 16mm.
- Balloon diameter – 16mm.
- Working length – 200mm.
- Compatible with 0,035” wire guide.
- Catheter diameter – 7Fr.
- Individually packed for **Single use only**.



Ordering information – Mention product code and quantity required

Reference	Catheter Diameter (mm/Fr)	Functional Catheter Length	Diameter of Inflated Balloon	Balloon Max Gas Installation	Balloon Max Fluid Injection	Recommended guiding wire	Guiding wire supplied	Syring CC
SSEC - 7	2.35/7	200mm	16mm	4.5cm ²	3.0ml	0.035”	NON	5

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises