


Wire Guided Monofilament Memory Basket

Used to retrieve foreign body/ stone from G.I tract – through the scope.

Technical information:

- Two basket sizes to facilitate stone retrieval from biliary tract.
- Basket shapes – Oval, Hexagonal.
- Retains shape for a longer time.
- Handle - light weight with flush port for contrast.
- Centre port for guide wire.
- **Sterile for Single use only.**



NOT RECOMMENDED FOR LITHOTRIPSY

Ordering information – Mention product code and quantity required

Reference	Product Name	Required Channel	Basket	Working Length
SMMB – WG - 30	Monofilament Memory Basket Wire Guided	2.8mm	30mm	230cm
SMMB – WG - 35	Monofilament Memory Basket Wire Guided	2.8mm	35mm	230cm
SMMB – WG - 40	Monofilament Memory Basket Wire Guided	2.8mm	40mm	230cm

Contact:-

www.makendoscopy.com

mak_enter@yahoo.com

mak_enter62@rediffmail.com

Phone - +91 22 25289181

MAK Enterprises