

## Biliary Dilator

**Used to dilate papilla or biliary strictures Through the Scope.**

### Technical Information

- Tapering length - 3 cm.
- Size - 5fr. - 6fr. - 7 fr.- 8.5 fr. - 9fr. - 10fr. - 11.5fr.
- Radioopaque PTFE material .
- Availability - Sterile - Single use only.
- Individually packed.



**Ordering Information**-Mention product code x quantity required

Reference	Product Name	TipØ	Required ChannelØ	Working Length
S - BD- 5	Biliary Dilator 5 Fr.	2.5 Fr.	3.2 mm	200 cm.
S - BD- 6	Biliary Dilator 6 Fr.	3.5 Fr.	3.2 mm	200 cm.
S - BD- 7	Biliary Dilator 7 Fr.	4 Fr.	3.2 mm	200 cm.
S - BD- 8.5	Biliary Dilator 8.5 Fr.	5.5 Fr.	3.2 mm	200 cm.
S - BD- 9	Biliary Dilator 9 Fr.	6 Fr.	3.2 mm	200 cm.
S - BD- 10	Biliary Dilator 10 Fr.	6 Fr.	4.2 mm	200 cm.
S - BD- 11.5	Biliary Dilator 11.5 Fr.	8.5 Fr.	4.2 mm	200 cm.
S - BD- 5-DBE	Biliary Dilator 5 Fr. For Double Balloon Enteroscope	2.5 Fr.	3.2 mm	280 cm.
S - BD- 6-DBE	Biliary Dilator 6 Fr. For Double Balloon Enteroscope	3.5 Fr.	3.2 mm	280 cm.
S - BD- 6-DBE	Biliary Dilator 6 Fr. For Double Balloon Enteroscope	3.5 Fr.	3.2 mm	280 cm.
S - BD- 7-DBE	Biliary Dilator 7 Fr. For Double Balloon Enteroscope	4 Fr.	3.2 mm	280 cm.
S - BD- 8.5-DBE	Biliary Dilator 8.5 Fr. For Double Balloon Enteroscope	5.5 Fr.	3.2 mm	280 cm.
S - BD- 9-DBE	Biliary Dilator 9 Fr. For Double Balloon Enteroscope	6 Fr.	3.2 mm	280 cm.
S - BD- 10-DBE	Biliary Dilator 10 Fr. For Double Balloon Enteroscope	6 Fr.	4.2 mm	280 cm.
S - BD-11.5-DBE	Biliary Dilator 11.5 Fr. For Double Balloon Enteroscope	8.5 Fr.	4.2 mm	280 cm.