

## Biopsy Forcep (Broncho)

Used endoscopically to obtain mucosal tissue samples from bronchiole for microscopic histopathological examination.

### Technical Information :

- Availability - Sterile For single use only.
- Compatible With - 2.0 / 2.8 mm instrument channel
- Length : 110 Cm
- Individually Packed



**Ordering Information-** Mention the product code x quantity required.

Reference	Product Name	Type	Required Channel Ø	Working length
S- BF- B	Biopsy Forcep	Without Spike	2.0/2.8 mm	110 cm