

Cytology Brush (Gastro)

Used to collect cells from the upper and lower Gastrointestinal Tracts

Technical Information

- Length - 170cm / 230cm & 325cm
- Wire Guide Compatibility - 0.035"
- Availability - Individually packed. Sterile For single use only.



Ordering Information-Mention product code x quantity required

Reference	Product Name	Required Channel Ø	Length
S - CY-170	CYTOLOGY BRUSH	2.8 mm	170 cm
S - CY-230	CYTOLOGY BRUSH	2.8 mm	230 cm
S - CY-DL	WIRE GUIDED CYTOLOGY BRUSH	230 cm	USH
S - CY - S	CYTOLOGY BRUSH	2.8 mm	1.2 mm