

Endoscopic Scissor

Used to cut and coagulate during endoscopic mucosal resection.

Technical Information:

- Working length 160/230cm,
- Sheath diameter 2.5 mm
- Availability - Sterile for single use only.
- **Rotatable scissor is available**
- Active Cord available separately



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Length	Sheath Diameter
S-ES-160-H	ENDO SCISSOR HOT	160cm	2.5mm
S-ES-160-R-H	ENDO SCISSOR HOT ROTATABLE	160cm	2.5mm
S-ES-230-H	ENDO SCISSOR HOT	230cm	2.5mm
S-ES-230-R-H	ENDO SCISSOR HOT ROTATABLE	230cm	2.5mm