

## Naso Biliary Drainage Catheter

**Used to drain bile from obstructed biliary duct through Naso.**

### Technical Information

- Size - 6fr / 7fr. / 8.5fr. /10fr.
- Working Length - 250cm
- Sterile for single use only.



**Ordering Information**-Mention product code x quantity required

Reference	Product Name	Size
S-NBDC-6	Naso Biliary Drainage catheter	6 fr
S-NBDC-7	Naso Biliary Drainage catheter	7 fr
S-NBDC-8.5	Naso Biliary Drainage catheter	8.5 fr
S-NBDC-10	Naso Biliary Drainage catheter	10 fr