


Needle Knife

Used for Precut Papillotomy.

Technical Information:

- Tapered tip.
- Individually packed.
- Supplied Sterile and it is for single use.
- Active Cord available separately.
- Wire Guided Needle Knife accept 0.035" Guide Wire



Ordering Information-Mention the product code x quantity required.

| Reference | Product Name | Length | Needle Length | Outer Ø |
|----------------|--|--------|---------------|---------|
| SNK-100 | Needle Knife | 200 cm | 5 mm | 1.7 mm |
| SNK-100-DL | Wire Guided Needle Knife | 200 cm | 5 mm | 2.45 mm |
| SNK-100-DBE | Needle Knife Double Balloon Endoscope | 208 cm | 5 mm | 1.7 mm |
| SNK-100-DL-DBE | Needle Knife Double Balloon Endoscope | 280cm | 5 mm | 2.45mm |