


## Neg Guide Wire (Uro)

Used to gain access Urology tract for access and exchange of accessories.

### Technical Information

- Size Diameter - 0.025 inch / 0.032 inch / 0.035 inch. / 0.038 inch.
- Availability - Individually packed.
- Sterile For single use only.
- Length - 150 cm.
- Tip - Straight / Angle / J



Ordering Information :

Reference	Product Name
S-NEG-025-S-150	Neg Guide Wire 0.025 Straight
S-NEG-035-S-150	Neg Guide Wire 0.035 Straight
S-NEG-025-A-150	Neg Guide Wire 0.025 Angulated
S-NEG-035-A-150	Neg Guide Wire 0.035 Angulated
S-NEG-025-J-150	Neg Guide Wire 0.025 J Tip
S-NEG-035-J-150	Neg Guide Wire 0.035 J Tip
S-NEG-032-S-150	Neg Guide Wire 0.032 Straight
S-NEG-038-S-150	Neg Guide Wire 0.038 Straight
S-NEG-032-A-150	Neg Guide Wire 0.032 Angulated
S-NEG-038-A-150	Neg Guide Wire 0.038 Angulated
S-NEG-032-J-150	Neg Guide Wire 0.032 J Tip
S-NEG-038-J-150	Neg Guide Wire 0.038 J Tip