


Single Action Stent Introducer Set

Used For Endoscopically Biliary Stent Placement to Reduce time and Easier Stent Placement.

Technical Information

- Size - 8.5fr., 10fr. & 11.5fr. Single action stent deployment.
- Radioopaque Marker Over Guiding Catheter.
- Low Friction Movement.
- Compatible With 0.035 Inch Ø Wire Guide.
- Pusher length - 180cm
- Guiding catheter length - 220cm
- Supply with Biliary Stent .
- Supply without Guide Wire.
- Availability - Sterile, For Single Use Only.
- Biliary Stent Length - 3 cm to 20 cm



Ordering Information- Mention the product code x required Length x QTY required.

Example : single action stent introducer system set X 10 fr X Teflon X 10 fr x Straight X 10cm long = SSISS-10-TEF-10-S-10

Specifications	Single Action Stent Introducer Set 8.5 Fr.	Single Action Stent Introducer Set 10 Fr.	Single Action Stent Introducer Set 11.5 Fr.
Material			
PTFE	SSISS-8.5-TEF-8.5	SSISS-10-TEF-10	SSISS-11.5-TEF-11.5
PE	SSISS-8.5-PU-8.5	SSISS-10-PU-10	SSISS-11.5-PU-11.5
Variant :			
1- Straight	SSISS-8.5-TEF-8.5-S SSISS-8.5-PU-8.5-S	SSISS-10-TEF-10-S SSISS-10-PU-10-S	SSISS-11.5-TEF-11.5-S SSISS-11.5-PU-11.5-S
2-Double Pig Tail	SSISS-8.5-PU-8.5-DP	SSISS-10-PU-10-DP	SSISS-11.5-PU-11.5-DP
3-Standard	SSISS-8.5-TEF-8.5-STD SSISS-8.5-PU-8.5-STD	SSISS-10-TEF-10-STD SSISS-10-PU-10-STD	SSISS-11.5-TEF-11.5-STD SSISS-11.5-PU-11.5-STD
4-Amsterdam	SSISS-8.5-TEF-8.5-AMST SSISS-8.5-PU-8.5-AMST	SSISS-10-TEF-10-AMST SSISS-10-PU-10-AMST	SSISS-11.5-TEF-11.5-AMST SSISS-11.5-PU-11.5-AMST
5- Pancreatic	SSISS-8.5-TEF-8.5-PAN SSISS-8.5-PU-8.5-PAN	SSISS-10-TEF-10-PAN SSISS-10-PU-10-PAN	SSISS-11.5-TEF-11.5-PAN SSISS-11.5-PU-11.5-PAN
6- S Banded	SSISS-8.5-PU-8.5-SB	SSISS-10-PU-10-SB	SSISS-11.5-PU-11.5-SB

7- Single Pigtail	SSISS-8.5-PU-8.5-SP	SSISS-10-PU-10 -SP	SSISS-11.5-PU-11.5 - SP
8- C Shape	SSISS-8.5-PU-8.5-CS	SSISS-10-PU-10 -CS	SSISS-11.5-PU-11.5 - CS