


Stent Pusher

Endoscopically biliary stent placement

Technical Information

- Size - 5 fr. - 7 fr. - 8.5 fr. - 10fr. - 11.5fr.
- Availability - Sterile - Single use only.
- Individually packed.
- Length as per customer requirement



Ordering Information-Mention the product code x quantity required.X Length as per customer requirement

Reference	Product Name	Size	Required ChannelØ
SSP-5 -X Length	Biliary Stent Pusher	5 Fr.	≥ 3.2 mm
SSP-7 -X Length	Biliary Stent Pusher	7 Fr.	≥ 3.2 mm
SSP-8.5 -X Length	Biliary Stent Pusher	8.5 Fr.	4.2 mm
SSP-10 -X Length	Biliary Stent Pusher	10 Fr.	4.2 mm
SSP-11.5 -X Length	Biliary Stent Pusher	11.5 Fr.	4.2 mm