


Stone Basket

Indicated for the manipulation and removal / retrieval of ureteral stones from the bladder, renal pelvis and ureter.

Technical Information

- Sheath Diameter - 1.8 fr, 2.4 fr, 3 fr, 4 fr & 5 fr
- Availability - Individually packed.
- Sterile For single use only.
- Length - Custom
- Type - Helical = H / Tipless = T / Straight = F / Tipless Helical = TH
- Wire Qty. - 3/4/5 wires.
- **With Dismantle facility**



Ordering Information :

Reference : SEMB - Size Diameter - Wire Qty. - Length - Shape

For Example : If you want 3fr., 4 wire, 90cm, Tipless

Order Ref. No. is SEMB - 3490 - T

Tailormade device can be produced provided by pre approved drawing/ specification