


Stone Extraction Balloon

Used For Endoscopic Extraction of Biliary Stones .

Technical Information:

- Triple lumen feature provides wire guide access to facilitate stone extraction, from desired location during procedure.
- Catheter Diameter-7 Fr .
- Compatible with 0.035 inch wire guide.
- Availability-Individually packed. Sterile for single use only.
- Working length-200 cm.
- Balloon diameter-16 mm.
- Balloon length-16 mm.



Ordering Information-Mention product code x quantity required

Reference	Catheter Diameter FR	Functional Catheter Length	Diameter of Inflated Balloon	Balloon Max Air Inflation	Balloon Max Fluid Injection	Recommended Guiding Wire(Inch)	Syringe CC	Tip Diameter FR
SSEC - 7	7	200 cm	16mm	3.5 cc	3.0 ml	0.035"	5	7