


Stone Grasper

Used to facilitate endoscopic manipulation / removal of stones and retrieval of foreign objects. Can also be used for removing stents through either rigid or flexible cystoscopy.

Technical Information

- Sheath Diameter - 3 fr / 4 fr / 5 fr / 7 fr
- Availability - Individually packed.
- Sterile For single use only.
- Length - As per Customer Requirement
- Type - 3 Prong



Ordering Information : Mention the product code x quantity required.

Reference	Product Name	Size
S-SG-3XLENGTH	Stone Grasper - 3 Nail	3 fr
S-SG-4XLENGTH	Stone Grasper - 3 Nail	4 fr
S-SG-5XLENGTH	Stone Grasper - 3 Nail	5 fr
S-SG-7XLENGTH	Stone Grasper - 3 Nail	7 fr

Tailormade device can be produced provided by pre approved drawing/ specification