


Aspiration Needle

Used for endoscopic aspiration from lungs and bronchiole

Technical Information -

- Sharp and short bevel. 19 / 21 swg
- Working length 110cm
- Compatible with 2.8 mm instrument channel
- Availability - sterile for single use only.
- Needle length 15 mm



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Needle Gauge	Working Length	Required Channel Ø
SAN - 19	ASPIRATION NEEDLE	19 G	110 cm	2.8 mm.
SAN - 21	ASPIRATION NEEDLE	21 G	110 cm	2.8 mm.