

## Biliary Balloon Dilator

Used to dilate papilla or biliary strictures through the scope

### Technical Information:

- Inflated balloon diameter - 6mm,8mm, 10 mm
- Balloon length : 4 cm
- Shaft diameter : 2.5 mm
- Working length : 2000 mm
- Wire guide compatibility : 0.035"
- Atraumatic Tip.
- Sterile for single use only.



Ordering Information-Mention product code x quantity required

Reference	Product Name	BalloonØ	Burst pressure
S-BBD-640	Biliary Balloon Dilator	6 mm	10ATM
S-BBD-840	Biliary Balloon Dilator	8 mm	10ATM
S-BBD-1040	Biliary Balloon Dilator	10 mm	10ATM