


Biliary Metal Stent

Used for Stenosis or Obstructed Biliary tract caused by benign or malignant changes.

Technical Information

- With Introducer System.
- Radiotherapy compatible.
- Expanded Stent Diameter. - 10mm.
- Stent Length - 60mm / 80mm / 100mm/120mm/150mm
- Uncovered stent.
- Compatible with 0.038 inch wire guide.
- Compatible with instrument channel size - 3.2mm
- Availability - Sterile - For single use only.



Ordering Information-Mention the product code x quantity required.

Reference	Product Name	Stent Length	Stent Ø
S - MBS - 6	Biliary Metal Stent	60mm	10mm
S - MBS - 8	Biliary Metal Stent	80mm	10mm
S - MBS - 10	Biliary Metal Stent	100mm	10mm
S - MBS - 12	Biliary Metal Stent	120mm	10mm
S - MBS - 15	Biliary Metal Stent	150mm	10mm