


Biliary Stent

Used to drain bile from obstructed biliary duct.

Technical Information

- Diameter 5 Fr. - 7Fr. - 8.5 Fr. - 10 Fr. & 11.5 Fr. Size
- Length 3 cm. to 20 cm. as per customer requirement.
- Material choice- Radioopaque PTFE OR PE.
- Varients
 - Straight-S
 - Double Pigtail-DP
 - Amsterdam-AMST
 - Standard-STD
 - Pancreatic-PAN
 - S Banded-SB
 - Single Pigtail-SP
 - C Shape-CS
 - Tree Type (Pine tree)-T
- Available-Sterile, Single use.
- Wireguide Compatibility-0.035'.



Ordering Information- Mention the product code x required Length x QTY required.

Example : Biliary Stent X Plastic X 7fr X Standard X10cm length = SBS-PU-7-STD-10

Specifications	Stent 5 Fr.Reference	Biliary Stent 7 Fr.Reference	Biliary Stent 8.5 Fr.Reference	Biliary stent 10Fr.Reference	Biliary stent 11.5 Fr.Reference
Material					
PTFE	SBS-TEF-5	SBS-TEF-7	SBS-TEF-8.5	SBS-TEF-10	SBS-TEF-11.5
PE	SBS-PU-5	SBS-PU-7	SBS-PU-8.5	SBS-PU-10	SBS-PU-11.5
Varient :					
1-Straight	SBS-TEF-5-S	SBS-TEF- 7- S	SBS-TEF-8.5-S	SBS-TEF-10-S	SBS-TEF-11.5-S
	SBS-PU-5-S	SBS-PU- 7-S	SBS-PU-8.5-S	SBS-PU-10-S	SBS-PU-11.5-S
2-Double Pig Tail	SBS-PU-5-DP	SBS-PU- 7-DP	SBS-PU-8.5-DP	SBS-PU-10-DP	SBS-PU-11.5-DP
3-Standard	SBS-TEF-5-STD	SBS-TEF-7-STD	SBS-PU-8.5-STD	SBS-TEF-10-STD	SBS-TEF-11.5-STD
	SBS-PU-5-STD	SBS-PU-7-STD	SBS-TEF-8.5-STD	SBS-PU-10-STD	SBS-TEF-11.5-STD
	SBS-TEF-5-AMST	SBS-TEF-7-AMST	SBS-TEF-8.5-AMST	SBS-TEF-10-AMST	SBS-TEF-11.5-AMST

4-Amsterdam	SBS-PU-5-AMST	SBS-PU-7-AMST	SBS-PU-8.5-AMST	SBS-PU-10-AMST	SBS-PU-11.5-AMST
5- Pancreatic	SBS-TEF-5-PAN	SBS-PU-7-PAN	SBS-TEF-8.5-PAN	SBS-PU-10-PAN	SBS-TEF-11.5-PAN
	SBS-PU-5-PAN	SBS-TEF-7-PAN	SBS-PU-8.5-PAN	SBS-TEF-10-PAN	SBS-PU-11.5-PAN
6- S Banded	SBS-PU-5-SB	SBS-PU-7-SB	SBS-PU-8.5-SB	SBS-PU-10-SB	SBS-PU-11.5-PAN
7- Single Pigtail	SBS-PU-5-SP	SBS-PU-7-SP	SBS-PU-8.5-SP	SBS-PU-10-SP	SBS-PU-11.5-SP
8- C Shape	SBS-PU-5-CS	SBS-PU-7-CS	SBS-PU-8.5-CS	SBS-PU-10-CS	SBS-PU-11.5-SP
	SBS-TEF-5-CS	SBS-TEF-7-CS	SBS-TEF-8.5-CS	SBS-TEF-10-CS	SBS-TEF-11.5-CS
9- Tree Type			SBS-TEF-8.5-T	SBS-TEF-10-T	SBS-TEF-11.5-T