


Cleaning Brush (Broncho)

Used to clean endoscope channel.

Technical Information

- Type both end bristles
- Shaft plastic coated wire
- Semi soft and channel friendly bristles for cleaning of endoscope channel.
- Atraumatic tip
- Compatible with 2.0 / 2.8 mm instrument channel
- Length 110cm
- Individually packed



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Required Channel Ø	Working Length
S-CB-BR	Cleaning Brush Single use	2.0 mm.	110 cm
S-CB-BR-R	Cleaning Brush Reusable	2.8 mm.	110 cm