

ERPC Cannula

Used endoscopically to opacify the biliary duct with radioopaque contrast Material for fluoroscopic view.

Technical Information

- Type - Tapered tip / Ball tip / Double lumen
- Catheter tip Diameter - 4.5 Fr. / 5.5 Fr.
- Recommended wire guide Diameter 0.035 inch.
- Availability - Sterile For single use only.
- Length- 200 cm & 280 cm
- Individually packed.



Ordering Information-Mention product code x quantity required

Reference	Product Name	TipØ	Working Length	Required ChannelØ
S- ERC-T-SL	Single Lumen ERCP Cannula Tapered Tip	4.5 Fr	200 cm	2.8 mm
S- ERC-B-SL	Single Lumen ERCP Cannula Ball Tip	5.5 Fr	200 cm	2.8 mm
S- ERC-T-DL	Single Lumen ERCP Cannula Tapered Tip	5.5 Fr	200 cm	2.8 mm
S- ERC-T-SL-DBE	Single Lumen ERCP Cannula Tapered Tip For Double Balloon Enteroscope	4.5 Fr	280 cm	2.8 mm
S- ERC-B-SL-DBE	Single Lumen ERCP Cannula Ball Tip For Double Balloon Enteroscope	5.5 Fr	280 cm	2.8 mm
S- ERC-T-DL-DBE	Double Lumen ERCP Cannula Tapered Tip For Double Balloon Enteroscope	5.5 Fr	280 cm	2.8 mm