

Esophageal Balloon Dilator

Used to dilate esophageal strictures over previously placed wire guide through the scope

Technical Information:

- Inflated balloon diameter - 6mm, 8mm, 10 mm, 12 mm, 14 mm, 16mm, 18mm
- Shaft diameter : 2.5 mm
- Working length : 1800 mm
- Wire guide compatibility 0.035"
- Sterile for single use only.



Ordering Information-Mention product code x quantity required

Reference	Product Name	Balloon Length	Inflated Balloon Ø	Burst Pressure
S-EBD-640	Esophageal Balloon Dilator	40mm	6mm	10 ATM
S-EBD-840	Esophageal Balloon Dilator	40 mm	8mm	10 ATM
S-EBD-1040	Esophageal Balloon Dilator	40 mm	10mm	10 ATM
S-EBD-1260	Esophageal Balloon Dilator	60 mm	12mm	6 ATM
S-EBD-1460	Esophageal Balloon Dilator	60 mm	14mm	6 ATM
S-EBD-1680	Esophageal Balloon Dilator	80 mm	16mm	6 ATM
S-EBD-1880	Esophageal Balloon Dilator	80 mm	18mm	5 ATM