

## Esophageal Dilator

**Used to Dilate Esophageal Strictures over previously placed wire guide.**

### Technical Information

- Available Sizes - 5mm to 20mm Outer Diameter.
- Other sizes can be produced on specific request.
- Soft, Flexible and Radioopaque Tip & Radioopaque Junction.
- Gradual Tapering.
- Low friction lumen surface.
- Dilator Size-Diameter-5 mm. to 20 mm.Length 100 cm.,70 cm.

### Availability-

1. Set of 6 Pcs. Even diameter./Odd diameter.
2. Set of 12 Pcs. & 16 Pcs.
3. Non Sterile & Reusable

Set includes 1x Shaili Flexible spring tip Guide Wire 250 cm long. Carrying case is available with order of set only.



**Ordering Information-**Mention product code x quantity required

Reference	Outer Ø	Length
SED-500	5mm	100cm
SED-600	6mm	100cm
SED-700	7mm	100cm
SED-800	8mm	100cm
SED-900	9mm	100cm
SED-1000	10mm	100cm
SED-1100	11mm	100cm
SED-1200	12mm	100cm
SED-1300	13mm	100cm
SED-1400	14mm	100cm
SED-1500	15mm	100cm

SED-1600	16mm	100cm
SED-1700	17mm	100cm
SED-1800	18mm	100cm
SED-1900	19mm	100cm
SED-2000	20mm	100cm

Reference	Outer Ø	Length
SED-570	5mm	70cm
SED-670	6mm	70cm
SED-770	7mm	70cm
SED-870	8mm	70cm
SED-970	9mm	80cm
SED-1070	10mm	70cm
SED-1170	11mm	70cm
SED-1270	12mm	70cm
SED-1370	13mm	70cm
SED-1470	14mm	70cm
SED-1570	15mm	70cm
SED-1670	16mm	70cm
SED-1770	17mm	70cm
SED-1870	18mm	70cm
SED-1970	19mm	70cm
SED-2070	20mm	70cm