

Foreign Body Holding Forcep 2 Prong, 3 Prong, Memory Basket (Broncho)

Used to retrieve foreign body from bronchiole tract.

Technical Information :

- Compatible with 2.0 / 2.8 mm instrument channel
- Working length 110cm
- Sterile for single use only



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Working Length	Required Channel Ø	Jaws
SFBH- 218-B	2 prong	110 cm.	2.0 mm	2 prong
SFBH- 318-B	3 prong	110 cm.	2.0 mm	3 prong
SMMB- BR	Memory Basket	110 cm.	2.0 mm	Basket