


## Hydra Neg Guide Wire (Uro)

Used to gain access Urology tract for access and exchange of accessories.

### Technical Information :

- Size Diameter - 0.025 inch. / 0.032 inch / 0.035 inch / 0.038 inch.
- Length - 150 cm.
- Tip - Straight / Angle / J, with Extra Filifom Wire
- Availability - Individually packed.
- Sterile For single use only.



### Ordering Information :

Reference	Product Name	Wire Ø
S-HYNEG-025-S-150	HYDRA NEG GUIDE WIRE	0.025 STRAIGHT
S-HYNEG-025-A-150	HYDRA NEG GUIDE WIRE	0.025 ANGULATED
S-HYNEG-025-J-150	HYDRA NEG GUIDE WIRE	0.025 J TIP
S-HYNG-032-S-150	HYDRA NEG GUIDE WIRE	0.032 STRAIGHT
S-HYNEG-032-A-150	HYDRA NEG GUIDE WIRE	0.032 ANGULATED
S-HYNEG-032-J-150	HYDRA NEG GUIDE WIRE	0.032 J TIP
S-HYNEG-035-S-150	HYDRA NEG GUIDE WIRE	0.035 STRAIGHT
S-HYNEG-035-A-150	HYDRA NEG GUIDE WIRE	0.035 ANGULATED
S-HYNEG-035-J-150	HYDRA NEG GUIDE WIRE	0.035 J TIP
S-HYNEG-038-S-150	HYDRANEG GUIDE WIRE	0.038 STRAIGHT
S-HYNEG-038-A-150	HYDRA NEG GUIDE WIRE	0.038 ANGULATED
S-HYNEG-038-J-150	HYDRA NEG GUIDE WIRE	0.038 J TIP