


Hyneg Combo Guide Wire (Uro)

Used to gain access Urology tract for access and exchange of accessories.

Technical Information

- Size Diameter - 0.032 inch / 0.035 inch.
- Availability - Individually packed.
- Sterile For single use only.
- One end with black hydrophilic tip / other end with exchange wire tip.
- Length - 150 cm
- Tip - Straight



Ordering Information : Mention the product code x quantity required.

Reference	Product Name
S-HYCG-032-S-150	Hyneg Combo Guide Wire 0.032 straight
S-HYCG-035-S-150	Hyneg Combo Guide Wire 0.035 straight