


## Single Action Stent Introducer System

**Used For Endoscopically Biliary Stent Placement to Reduce time and Easier Stent Placement.**

### Technical Information

- Size - 8.5fr., 10fr. & 11.5fr. Single action stent deployment.
- Radioopaque Marker Over Guiding Catheter.
- Low Friction Movement.
- Compatible With 0.035 Inch Ø Wire Guide.
- Supply Without Guide Wire and Stent .
- Availability - Sterile, For Single Use Only.
- Pusher length - 180cm
- Guiding catheter length - 220cm



Ordering Information- Mention the product code x required Length x QTY required.

Reference	Product Name	Size	Remark
SSIS - 8.5	Single Action Stent Introducer System	8.5 Fr.	Without Guide Wire and Stent.
SSIS - 10	Single Action Stent Introducer System	10 Fr.	Without Guide Wire and Stent.
SSIS - 11.5	Single Action Stent Introducer System	11.5 Fr.	Without Guide Wire and Stent.