


## Sphincterotome

**Used for Endoscopic Cannulation of Biliary Ductal System and Sphincterotomy.**

### Technical Information:

- Triple Lumen with Precurved tapered tip.
- Tip diameter 5.5fr.
- Working length 200 cm.
- Wire Guide Compatibility 0.035"
- Cutting wire with or without insulation as required
- Protective insulation sheath, unique protective feature for endoscope,
- Insulated wire length - for 20mm Loop - 4 mm, 25mm Loop - 8mm, 30 mm Loop - 10mm
- Supplied Sterile for single use.
- Active Cord available separately.



**Ordering Information**-Mention the product code x quantity required.

Reference	Product Name	Tip Length	Total Loop	With/Without Insulation
SSPH-SN-20	Sphincterotome Short Nose	5.0mm	20mm	With Insulation
SSPH-SN-25	Sphincterotome Short Nose	5.0mm	25mm	With Insulation
SSPH-SN-30	Sphincterotome Short Nose	5.0mm	30mm	With Insulation
SSPH-SN-NI-20	Sphincterotome Short Nose	5.0mm	20mm	Without Insulation
SSPH-SN-NI-25	Sphincterotome Short Nose	5.0mm	25mm	Without Insulation
SSPH-SN-NI-30	Sphincterotome Short Nose	5.0mm	30mm	Without Insulation
SSPH-LN-20	Sphincterotome Long Nose	8.0mm	20mm	With Insulation
SSPH-LN-25	Sphincterotome Long Nose	8.0mm	25mm	With Insulation
SSPH-LN-30	Sphincterotome Long Nose	8.0mm	30mm	With Insulation
SSPH-LN-NI-20	Sphincterotome Long Nose	8.0mm	20mm	Without Insulation
SSPH-LN-NI-25	Sphincterotome Long Nose	8.0mm	25mm	Without Insulation
SSPH-LN-NI-30	Sphincterotome Long Nose	8.0mm	30mm	Without Insulation
SSPH-NN-20	Sphincterotome No Nose	0.0mm	20mm	With Insulation
SSPH-NN-25	Sphincterotome No Nose	0.0mm	25mm	With Insulation

SSPH-NN-30	Sphincterotome No Nose	0.0mm	30mm	With Insulation
SSPH-NN-NI-20	Sphincterotome No Nose	0.0mm	20mm	With Insulation
SSPH-NN-NI-25	Sphincterotome No Nose	0.0mm	25mm	With Insulation
SSPH-NN-NI-30	Sphincterotome No Nose	0.0mm	30mm	With Insulation
SSPH-RC-NI-20	Sphincterotome Reverse Cut	2.0mm	20mm	Without Insulation
SSPH-RC-NI-25	Sphincterotome Reverse Cut	2.0mm	25mm	Without Insulation
SSPH-RC-NI-30	Sphincterotome Reverse Cut	2.0mm	30mm	Without Insulation
SSPH-SN-20-DBE	Sphincterotome Short Nose 280 cm.long	5.0mm	20mm	With Insulation
SSPH-SN-25-DBE	Sphincterotome Short Nose 280 cm.long	5.0mm	25mm	With Insulation
SSPH-SN-30-DBE	Sphincterotome Short Nose 280 cm.long	5.0mm	30mm	With Insulation
SSPH-LN-20-DBE	Sphincterotome Long Nose 280 cm.long	8.0mm	20mm	With Insulation
SSPH-LN-25-DBE	Sphincterotome Long Nose 280 cm.long	8.0mm	25mm	With Insulation
SSPH-LN-30-DBE	Sphincterotome Long Nose 280 cm.long	8.0mm	30mm	With Insulation