

Tip Guided Retrieval Basket


Used to retrieve Foreign Bodies and C.B.D. stones endoscopically

Technical Information:

- Working length 230 cm.
- Wire access point on tip of basket.
- Compatible with instrument channel 2.8 mm or more.
- Basket diameter -15mm. / 20 mm. / 25mm. / 30 mm. / 35 mm. / 40 mm. / 45 mm. / 50 mm. Multifilament, Hexagonal shaped.
- Flush port for contrast.

Compatible with Lithotripter.

- Sterile for single use only.



Ordering Information-Strand diameter is 0.41mm, Working length 230 cm.

Reference	Product Name	Basket Ø	Sheath Ø	Wire Guide Compitibility
SRB-WG-4152-T	Tip Guided Retrieval Basket	15 mm	2.5 mm	0.035"
SRB-WG-4202-T	Tip Guided Retrieval Basket	20mm	2.5 mm	0.035"
SRB-WG-4252-T	Tip Guided Retrieval Basket	25 mm	2.5 mm	0.035"
SRB-WG-4302-T	Tip Guided Retrieval Basket	30 mm	2.5 mm	0.035"
SRB-WG-4352-T	Tip Guided Retrieval Basket	35 mm	2.5 mm	0.035"
SRB-WG-4402-T	Tip Guided Retrieval Basket	40 mm	2.5 mm	0.035"
SRB-WG-4452-T	Tip Guided Retrieval Basket	45 mm	2.5 mm	0.035"
SRB-WG-4502-T	Tip Guided Retrieval Basket	50 mm	2.5 mm	0.035"

Tailormade device can be produced provided by pre approved drawing/ specification